

<https://github.com/luketschepen/EE2361Final/tree/main>

INTRO- Description of the Library and its overall functions

Members: Dhayalan Balasubramanian, Luke Tschepen, Nicholas Kang, Cole Jaeger.

The Blood Oxygen Sensor project integrates the MAX30102 sensor for pulse oximetry and heart-rate monitoring, alongside the Grove RGB LCD module for versatile display feedback. Physical input is managed by the `button_switch.c` file, handling interaction with a button switch for toggling display and capturing button events. The `algorithm.c` file computes heart rate and SpO2 levels using peak detection techniques and precise algorithms on sensor data. The `I2C.h` file handles essential I2C communication for interfacing with the sensor and other peripherals.

Separate Descriptions

MAX 30102 Description:

The code implements functionalities for interfacing with the MAX30102 sensor, which is a pulse oximeter and heart-rate sensor. It includes functions for initializing the sensor, configuring its settings, reading sensor data, processing the signals, and displaying the results.

The Grove RGB LCD module is a multifunctional display that integrates an RGB backlight and a character LCD. It is designed to provide a versatile display solution for various electronic projects. This document outlines the functionality and usage of the Grove RGB LCD module within the context of the Blood Oxygen Sensor project.

The `button_switch.c` file provides functionality for interfacing with a physical button switch in the Blood Oxygen Sensor project. It includes functions for initializing the button, detecting button presses, and determining the state of the button. The button switch is connected to RB7 (Pin 6) of the microcontroller, and its usage involves capturing button press events and toggling the display of heart rate information.

The `algorithm.c` file is a crucial component of the Blood Oxygen Sensor project, dedicated to computing heart rate and SpO2 levels. It implements peak detection techniques on data from IR and red sensors to identify PPG signal peaks for heart rate determination. Additionally, it employs a precise algorithm to derive SpO2 values from the ratio of AC to DC components in the IR and red signals.

The I2C.h file contains functions for initializing and managing I2C communication on PIC24 microcontrollers. It includes functions for initializing I2C parameters, initiating start and stop conditions, writing and reading data over I2C, handling acknowledgments, and performing repeated start conditions.

Hardware description - what device(s), part numbers, links, etc.

SNAP: MP LAB SNAP PG164100

<https://ww1.microchip.com/downloads/aemDocuments/documents/DEV/ProductDocuments/UserGuides/50002787C.pdf>

PIC24: PIC24FJ64GA002

<https://ww1.microchip.com/downloads/en/DeviceDoc/39881e.pdf>

Sensor: MAX30102

<https://www.analog.com/media/en/technical-documentation/data-sheets/max30102.pdf>

RGB LCD: Grove-LCD RGB Backlight JHD1313

https://files.seeedstudio.com/wiki/Grove_LCD_RGB_Backlight/res/JHD1313%20FP-RGB-1%201.4.pdf

Button: Square Tactile Pushbutton Switch SPST Momentary -ELSW-TPB-SPST

<https://www.ckswitches.com/media/1462/pts125.pdf>

Full documentation- description - what device(s), part numbers, links, etc.

MAX 30102 Description:

The code implements functionalities for interfacing with the MAX30102 sensor, which is a pulse oximeter and heart-rate sensor. It includes functions for initializing the sensor, configuring its settings, reading sensor data, processing the signals, and displaying the results.

Devices and Parts:

SNAP

PIC24

MAX30102 Sensor

RGB LCD display

Button switch

Links to datasheets and part numbers:

SNAP: MP LAB SNAP PG164100

<https://ww1.microchip.com/downloads/aemDocuments/documents/DEV/Product Documents/UserGuides/50002787C.pdf>

PIC24: PIC24FJ64GA002

<https://ww1.microchip.com/downloads/en/DeviceDoc/39881e.pdf>

Sensor: MAX30102

<https://www.analog.com/media/en/technical-documentation/data-sheets/max30102.pdf>

RGB LCD: Grove-LCD RGB Backlight JHD1313

https://files.seeedstudio.com/wiki/Grove_LCD_RGB_Backlight/res/JHD1313%20FP-RGB-1%201.4.pdf

Button: Square Tactile Pushbutton Switch SPST Momentary -ELSW-TPB-SPST

<https://www.ckswitches.com/media/1462/pts125.pdf>

Public Functions:

Initialization and Configuration:

max30102_init(void): Initializes the MAX30102 sensor.

max30102_setup_spo2(): Configures the MAX30102 sensor for SpO2 measurement.

max30102_setup_default(void): Sets up the MAX30102 sensor to the default state.

max30102_setup_display(void): Configures the display for showing SpO2 and heart rate data.

uint32_t max30102_read_partID(); Reads the Part ID(0x15) and was written to ensure the hardware, sensor, and I2C protocol were working.

Data Acquisition and Processing:

max30102_process_signals(): Processes the sensor signals to calculate SpO2 and heart rate.

max30102_readTemp(): Reads the die temperature from the sensor.

calculateMovingAverage(uint32_t *buffer, int size): Calculates the moving average of the data buffer.

Buffer Management:

max30102_RED_putBuff(long int data): Inserts data into the Red LED buffer.

max30102_IR_putBuff(long int data): Inserts data into the Infrared LED buffer.

max30102_RED_getBuff(): Retrieves data from the Red LED buffer.

max30102_IR_getBuff(): Retrieves data from the Infrared LED buffer.

Control and Configuration:

Functions to configure various parameters such as LED mode, ADC range, sample rate, pulse width, FIFO settings, etc.

Utility:

max30102_shutdown(void): Shuts down the sensor.

softReset(void): Performs a soft reset of the sensor.

Arguments and Outputs:

Many functions take parameters to configure specific settings on the sensor, such as LED mode, sample rate, pulse width, etc.

The functions for inserting and retrieving data from the buffers

(max30102_RED_putBuff, max30102_IR_putBuff, max30102_RED_getBuff, max30102_IR_getBuff) accept and return data values.

Functions for processing signals (max30102_process_signals()) calculate SpO2 and heart rate and store the results in global variables.

Display functions (max30102_setup_default, max30102_setup_display) output data to an RGB LCD display.

Grove RGB LCD Documentation

Description:

The Grove RGB LCD module is a multifunctional display that integrates an RGB backlight and a character LCD. It is designed to provide a versatile display solution for various electronic projects. This document outlines the functionality and usage of the Grove RGB LCD module within the context of the Blood Oxygen Sensor project.

Devices and Parts:

SNAP

PIC24 microcontroller

MAX30102 Sensor

RGB LCD display

Button switch

Links:

Public Functions:

// Function prototypes

```
void delay_ms(unsigned int ms); // Delays execution for a specified number of
milliseconds.
void rgb_cmd(char Package); // Sends a command to control RGB LEDs.
void setReg(unsigned char reg, unsigned char dat); // Sets a register with specified
data.
void setRGB(unsigned char r, unsigned char g, unsigned char b); // Sets RGB LED
color with specified red, green, and blue values.
void rgb_clr(); // Clears RGB LED color.
void rgb_home(); // Sets RGB LEDs to default state.
void noDisplay(); // Turns off display.
void display(); // Turns on display.
void blinkLED(void); // Blinks an LED.
void noBlinkLED(); // Stops LED blinking.
void setColorAll(); // Sets color for all LEDs.
void setColorWhite(); // Sets color to white for all LEDs.
void setColorRed(); // Sets color to red for all LEDs.
void setCursor(uint8_t x, uint8_t y); // Sets cursor position on a display.
void grovergb_init(void); // Initializes Grove RGB module.
void printStr(const char s[]); // Prints a string on a display.
void printChar(char myChar); // Prints a character on a display.
```

Initialization and Configuration:

grovergb_init(): Initializes the Grove RGB LCD module to be compatible with the PIC24 microcontroller using I2C communication.

Data Acquisition and Processing:

Not applicable in this document. (Related to the blood oxygen sensor functionality)

Buffer Management:

Not applicable in this document.

Control and Configuration:

setColorAll(): Turns off all colors on the display.
setColorWhite(): Sets the screen color to white.
setColorRed(): Sets the screen color to red.
blinkLED(): Flashes the specified LED registers.
noBlinkLED(): Stops the flashing of the LED registers.

Utility:

printStr(const char s[]): Prints a string to the display.
printChar(char myChar): Prints a character to the display.

Arguments and Outputs:

grovergb_init(): No arguments or outputs.

setColorAll(), setColorWhite(), setColorRed(), blinkLED(), noBlinkLED(): No arguments or outputs.

printStr(const char s[]): Takes a string s as input and prints it on the display.

printChar(char myChar): Takes a character myChar as input and prints it on the display.

Basic usage example- bare minimum to test hardware functionality

In order to test the functionality of the hardware the max30102_read_partID() function can be called. After calling this function the sensors will be configured and running and will display the part ID on the LCD.

You can also display the raw readings from the sensor using the following configuration.

```
int main(void) {
    setup(); // Perform all initial setup tasks

    while(1) {
        rgb_clr();
        getRead();
        printf("IR Reading: %u\n ", max30102_IR_getBuff());
        printf("RED Reading: %u\n", max30102_RED_getBuff());

    }
    return 0;
}
```

Advanced usage example - covering all the functions and feature

An advanced usage example could be displaying either heart rate or SpO2 and controlling the display with the press of a button similar to what we demo'd.

```
void max30102_process_signals(){
    getRead();
    // Read 100 samples into the buffers
    for (int i = 0; i < 100; i++) {
        irBuffer[i] = max30102_IR_getBuff(); // Get one IR sample
        redBuffer[i] = max30102_RED_getBuff(); // Get one Red sample
    }
}
```

```

// Process the collected data
maxim_heart_rate_and_oxygen_saturation(irBuffer, 100, redBuffer, &spo2, &validSpo2,
&heartRate, &validHeartRate);
if (spo2 > 85) {
    // Update heart rate buffer with latest value
    for (int i = 0; i < MOVING_AVG_SIZE_HR - 1; i++) {
        heartRateBuffer[i] = heartRateBuffer[i + 1];
    }
    if (heartRate > 0){
        heartRateBuffer[MOVING_AVG_SIZE_HR - 1] = heartRate;

    }

    // Update SpO2 buffer with latest value
    for (int i = 0; i < MOVING_AVG_SIZE_SPO2 - 1; i++) {
        spo2Buffer[i] = spo2Buffer[i + 1];
    }
    if (spo2 > 0){
        spo2Buffer[MOVING_AVG_SIZE_SPO2 - 1] = spo2;

    }

    // Calculate moving average of heart rate
    heartRateAvg = calculateMovingAverage(heartRateBuffer, MOVING_AVG_SIZE_HR);
    spo2Avg = calculateMovingAverage(spo2Buffer, MOVING_AVG_SIZE_SPO2);

    max30102_setup_display();

}
else {

    max30102_setup_default();
}
}

```

```

int main(void) {
    setup(); // Perform all initial setup tasks

```

```

while(1) {
    rgb_clr(); //clear screen
    max30102_process_signals(); //process readings and update display

}
return 0;
}

```

If you wanted to build upon this even further you could add the following. An even more advanced usage example could be implementing a power saving mode, By detecting how long it has been since the previous valid reading. If it has been over two minutes then go into “power saving mode” and the sensor would stop reading. The button could then be used to turn the sensor back on and continue reading. Another possible thing to add would be to specify your age and then corresponding “zones” would be set. Then based on the readings for heart rate and SpO2, we’d display either a red yellow or green backlight for how “healthy” the values are.

```

int main(void) {
    setup(); // Perform all initial setup tasks
    readForAge();//function to take in age and set zones. Uses button to take age.
    turnOnPowerSaving();//function to turn on power saving mode

while(1) {
    powerSaveCheck();//check for readings and if sensor should stay on/off
    rgb_clr();
    max30102_process_signals()
    }
return 0;
}

```


}

References used:

Basic Background Video - <https://www.youtube.com/watch?v=V5UvNVQsUsY>

Sensor:

Guide to using Sensor -

<https://www.instructables.com/Guide-to-Using-MAX30102-Heart-Rate-and-Oxygen-Sens/>

OxSensor Example Library

https://github.com/sparkfun/SparkFun_MAX3010x_Sensor_Library/blob/master/src/MAX30105.h

Oxsensor Arduino Tutorial

<https://lastminuteengineers.com/max30102-pulse-oximeter-heart-rate-sensor-arduino-tutorial/>

Data Sheet

<https://www.analog.com/media/en/technical-documentation/data-sheets/max30102.pdf>.

LCD:

Grove LCD wiki

https://wiki.seeedstudio.com/Grove-LCD_RGB_Backlight/

Grove LCD datasheet

https://files.seeedstudio.com/wiki/Grove_LCD_RGB_Backlight/res/JHD1313%20FP-RGB-1%201.4.pdf

Other LCD datasheet

https://files.seeedstudio.com/wiki/Grove-LCD_RGB_Backlight/Grove-LCD_RGB_Backlight_V5.0_Datasheet.pdf

Other Other LCD datasheet

https://files.seeedstudio.com/wiki/Grove_LCD_RGB_Backlight/res/PCA9633.pdf

PIC 24 Data Sheet:

<https://ww1.microchip.com/downloads/en/DeviceDoc/39881e.pdf>

PIC 24 I2C:

<https://ww1.microchip.com/downloads/en/devicedoc/70000195f.pdf>